Durevole

Pino de fibra de vidro com broca

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Fiber Glass Post with largo bur

USE INSTRUCTIONS

DESCRIPTION: DUREVOLE fiberglass posts are translucent, conical double intraradicular posts. The post has a modulus of elasticity similar to that of dentin, providing a more homogeneous mechanical unit, enabling a more uniform distribution of stress along the root. This package contains a drill for root canal preparation. For information on the measurements of the post and drill indicated, see the table 1.

RECOMMENDATION: Intraradicular posts are indicated to aid in the retention of extensive coronal restoration in endodontically treated teeth..

HOW TO USE:

1. Perform the radiographic examination (periapical radiography) and check the quality of the endodontic treatment and periapical health;
2. Select the post according to the diameter of the conduit positioning it on the radiograph of the tooth;
3. Perform absolute isolation of the operational field;
4. Remove endodontic filling material from 2/3 of the root canal using suitable instruments. Leave at least 3 to 5 mm of filling at the apex, ensuring apical sealing and avoiding bacterial contamination;
5. Select the width drill according to the post diameter, as indicated in the table at the end of the package insert. Prepare the root canal with wide drills at low speed and under refrigeration. The drill must be inserted into the root canal in the direction of the long axis of the canal, until it reaches, at most, the entire working part, according to the length to be prepared. Subsequently, the drill must be pressed laterally against the root dentin wall, along the entire contour of the canal. Avoid unnecessary wear;
6. Take a periapical radiograph of the selected post within the root canal. Check the image for the preparation of the conduit, the adaptation of the post and the condition of the remaining filling. The post must be perfectly adapted to the preparation, no spaces between the end of the post and the filling of the apical remnant;
7. If necessary, cut the coronal excess of the post with a diamond bur at high speed under refrigeration. Return post to position for fit check;
8. Clean the post with 70% alcohol and dry;
9. Apply a layer of silane to the entire surface of the post and allow it to dry for one minute. The post must be completely dry at the time of cementing;
10. Apply 37% phosphoric acid gel for 15 seconds into the root canal and coronary remnant. Wash with plenty of water and completely remove excess water from the canal with absorbent paper cones;
11. Apply the adhesive system (according to the manufacturer's instructions) along the entire length of the conditioned root canal and light cure;
12. Bring the resin cement into the root canal using appropriate instruments. If you prefer, cover the post with the luting agent as well;
13. Insert the post into the channel. Note that the post must be inserted to the bottom of the prepared cavity and the post must be completely covered by the resin cement;

14. Remove excess cement;
 15. Light cure the cement (according to the manufacturer's specification);
 16. It is recommended to perform a periapical radiograph to assess the final quality after cementation.
- COMPOSITION:** Posts: Fiberglass, epoxy resin and radiopacifying agent. Wide drill bits: stainless steel.
- POST CLEANING:** The cleaning of the post must be carried out with 70% alcohol.

BUR STERILIZATION: After using the bur, keep it immersed in enzymatic detergent. Subsequently, perform mechanical cleaning by friction using a brush and soap/detergent. Rinse thoroughly and dry the bur. Pack and/or condition the bur in material suitable for autoclave sterilization. The bur must be sterilized in an autoclave before use, at a temperature between 121°C and 135°C, observing the recommendations in the autoclave manufacturer's instruction manual. Repeat the process for each use.

PRECAUTIONS: The indication for the use of the post depends on a careful evaluation of the coronary remnant, the quality of the endodontic treatment and the periapical health of the tooth. The proper selection of drill diameter and post size is essential to avoid root fragility as well as to obtain the expected clinical results. Determine the proper working length and post length. Small ends of the fiberglass may be exposed and damage the skin. Post cutting can also expose sharp fibers. The bur must be used without excessive pressure on the dentinal walls. The sharp end of the drill bit can cause damage. The use of the drill, at low speed, and the cutting of the post with a diamond tip, at high speed, must be carried out under refrigeration. The use of dentin adhesives over the fiberglass post is unnecessary and not recommended. Use PPE to handle the product.

WARNINGS: This product is exclusively intended for dental use and must be handled in accordance with the instructions for use. The manufacturer is not responsible for damage caused by misuse or incorrect handling. Restricted use for dental professionals. The manipulation of the post must be done under absolute isolation. The fiberglass post should not be used in patients with hypersensitivity to the components of the product. The post is for single use and its reprocessing is prohibited. Over time, the wearability of the drill bit can be reduced. Do not use if there is loss of cutting ability. Do not use the product if it is out of date. For the disposal of the product, follow the legislation of your country. Keep out of the reach of children.


CONTRAINDICATIONS: The post is contraindicated for patients with parafunctional habits, patients with a coronary remnant smaller than 2 mm, patients whose endodontic treatment is unsatisfactory and/or patients with hypersensitivity to the components of the product formula. The drill should not be used on teeth with very thin dentin walls.

SHELF LIFE: 05 (five) years after the date of manufacture.

STORAGE: Store the product in its original, closed packaging, in a clean and dry place, with relative humidity below 95% and at a temperature of up to 50°C.

TABLE 1: MODELS

DUREVOLE MODEL N.	Colors	A Height	B Cervical Diameter	C Apical Diameter	BUR (N.)
0.3	Green	19mm	1,25	0,45	N/A
0.5	Purple	19mm	1,40	0,65	1
1.0	White	19mm	1,60	0,85	2
2.0	Yellow	19mm	1,80	1,05	3
3.0	Blue	19mm	2,00	1,25	4

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