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FDA Food and Drug Administration (USA).

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Endo ClarkTec Digital X-Ray Positioner

INSTRUCTIONS FOR USE

DESCRIPTION: Recommended for performing periapical radiographs of the upper and lower teeth (orthoradial, mesioradial and distoradial) for the purpose of radiographic diagnosis and intraoperative evaluation in endodontic treatment with the use of a digital sensor. Product designed in partnership with dental surgeon Andre Caliman.

COMPOSITION: Polyacetal.

INDICATION: Indicated to assist in the execution of periapical radiographs with the Clark technique (mesialization and distalization of the X-ray beam).

HOW TO USE:

To obtain x-rays during endodontic treatment of posterior teeth:

1. Select the radiographic positioner according to the tooth group to be x-rayed.
2. Assemble the radiographic positioner by fitting the ClarkTec device between the positioner body and the guidance ring and fitting the bite device next to the image receptor bracket.
3. Adjust the ClarkTec device according to the desired incidence (orthoradial, mesioradial or dystoradial). – See usage illustrations below.
4. Select the radiographic sensor holder according to your digital radiographic system and attach the holder to the digital sensor.
5. Attach the sensor-clamp assembly to the image receptor support.
6. Take the positioner to the patient's mouth, positioning the digital sensor so that it is centered on the tooth group to be x-rayed.
7. Ask the patient to bite down on the biting device to keep the positioner stable.
8. Align the cylinder of the x-ray machine with the orientation ring of the radiographic positioner.

To obtain x-rays during endodontic treatment of anterior teeth:

1. Select the radiographic positioner without the bite device.
2. Mount the radiographic positioner by fitting the ClarkTec device between the shaft and the guidance ring.
3. Adjust the ClarkTec device according to the desired incidence (orthoradial, mesioradial or dystoradial). – See usage illustrations below.
4. Select the radiographic sensor holder according to your digital radiographic system and attach the holder to the digital sensor.

5. Attach the sensor-clamp assembly to the image receptor support.
6. Take the positioner to the patient's mouth, positioning the digital sensor so that it is centered on the tooth group to be x-rayed.
7. Ask the patient to keep their mouth open and stabilize the positioner using the rod with the support of their hand.
8. Align the cylinder of the x-ray machine with the orientation ring of the radiographic positioner.

Note: The kit contains 18 fixing devices with different dimensions to adapt to the different digital sensors available on the market. Select the clamping device that is compatible with your digital sensor.

CLEANING: Use enzymatic detergent in ultrasonic tanks, according to manufacturers' instructions.

DISINFECTION: keep the radiographic positioner immersed in a 0.2% peracetic acid solution for 10 minutes. Perform visual inspection, dry the product with a clean cloth.

STERILIZATION: Pack the positioner in packaging suitable for sterilization in an autoclave at a temperature between 121°C and 135°C, as recommended by the autoclave manufacturer. The cleaning, disinfection and sterilization process must be carried out after each use.

DISPOSAL: The product must be disposed of in accordance with current local environmental legislation.

WARNINGS: Handle the device with care when moving the axis to the right or left.

PRECAUTIONS: Product indicated for professional use. Sterilize before use. Do not touch the product to the internal walls of the autoclave.

STORAGE: Store the product in its original, closed packaging, in a clean, dry place and away from direct light.

EXPIRY DATE: 05 years from the date of manufacture, in undamaged packaging. The useful life of this product depends on the care taken when handling it.

THIS PACKAGE CONTAINS:


1 piece "A" Endo RH Plus Digital; 1 piece "B" Endo RH Plus Digital; 2 orientation rings; 2 14 mm bite devices; 2 16 mm bite devices; 1 autoclavable pot; 18 digital radiographic sensor holders; 1 ClarkTec device.

ILLUSTRATIONS OF USE:

ASSEMBLY OF THE CLARKTEC DEVICE:


STEP 1:

Fit the ClarkTec device (2) next to the positioner body (3), in the position indicated below.


STEP 2:

Attach the orientation ring (1) to the other end of the ClarkTec device (2).



USE IN ORTHORADIAL INCIDENCE

Arrows on the device are aligned.

USE IN MESIORADIAL OR DYSTORADIAL INCIDENCE

The mesioradial or dystoradial incidence depends on the side of the dental arch.

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Rev 00 - 23.10.2023 | Cód 321510880012
ANVISA: 80213420034

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Medical device (Europe).

European Authorized Representative.

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