Sentinella – Indicador Biológico

IB.PHST08H – Leitura por pH – Vapor – 08 Horas

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IB.PHST08H – pH Indicator – Steam – 08 Hours

INSTRUCTIONS FOR USE:

DESCRIPTION: Self-contained biological indicator for steam sterilization with results in 08 hours.

INDICATION: Recommended for hospitals, clinics and healthcare institutions. Developed to monitor and ensure the effectiveness of vacuum-assisted sterilization and gravity air displacement at 121 – 135°C.

CONTENTS: Each vial contains a paper carrier with *Geobacillus stearothermophilus* spores packaged in a thermoplastic tube, along with a glass ampoule containing a pH-sensitive culture medium (purple color).

INSTRUCTIONS FOR USE:

1. Identify the biological indicator label with the autoclave identification number, load, and process date. Place it in an appropriate challenge device (DDP). Use another ampoule as a positive control in all tests, which should belong to the same batch and manufacturing date as the indicators subjected to the sterilization process. Identify the vial label with "C" (for "control") and the date.
2. Place the DDP along with the load in critical locations as specified by the autoclave manufacturer.
3. Sterilize according to the autoclave manufacturer's instructions.
4. After sterilization, wait for 10 minutes and remove the biological indicator from the device. **CAUTION:** The 10-minute time is important for the indicators to cool down and depressurize. Do not handle biological indicators without using safety equipment.
5. Check if the chemical process indicator printed on the biological indicator vial label has changed color (from pink to green), confirming exposure to the sterilization process. **IMPORTANT:** Color change is not evidence of a successful sterilization process. If the color has not changed, check the sterilization process.
6. Activate the biological indicator by breaking the glass ampoule, being careful not to break the plastic part. Check if the culture medium has been completely released onto the spore carrier.
7. Tap the bottom of the vial lightly with your finger to ensure the culture medium reaches the bottom of the vial, covering the paper carrier with the spores completely. **CAUTION:** Do not shake the ampoule up and down. The upper part has a filter that should not get wet.
8. Repeat the operation (starting from item 6) with the control vial that was not autoclaved.
9. Within two hours after sterilization, incubate the vials (test and control) at the appropriate temperature for the time specified on the product packaging.
10. The control ampoule should show a positive result (yellow color), and the test ampoule should show a negative result (purple color). The processed biological indicator results are not valid until the control ampoule reading is positive.

RESULTS ANALYSIS

Test	Control	Result
PURPLE	YELLOW	APPROVED
Expected Result (Positive): Indicates that there was no bacterial growth in the test ampoule and that there was bacterial growth in the control ampoule. Sterilization and testing were effective.		

Test	Control	Result
YELLOW	YELLOW	DISAPPROVED
Ambiguous Result (Negative): Indicates that sterilization was not effective – Sterilization process failure (check parameters and autoclave maintenance).		

Test	Control	Result
PURPLE	PURPLE	DISAPPROVED
Ambiguous Result (Negative): Indicates that there was no bacterial growth because the control ampoule was not activated, due to incorrect use or storage.		

Test	Control	Result
YELLOW	PURPLE	DISAPPROVED
Indicates a failure in the sterilization process and also that the control ampoule was not activated, indicating an error in use or incorrect storage.		

IMPORTANT: If a negative result is observed, repeat the test with a second biological indicator from the same batch as the first negative result.

STORAGE: Store at a temperature between 15°C and 30°C and relative humidity between 30 and 80%. Do not store the product near sterilizing agents such as ethylene oxide, dry heat, radiation, or other sterilization processes.

TRANSPORTATION: Transport in closed and reinforced boxes to avoid damage.

PRECAUTIONS: Incubation at 60°C should be performed in incubators compatible with the biological indicator. Single-use product. Make sure the ampoules were incubated for the specified time and temperature. If confirmed, the incubator may need maintenance. Incorrect storage conditions can also adversely affect biological indicators, even within their expiration date.

WARNINGS: If the sterilization process fails or does not reach the necessary sterilization parameters, the culture medium will change color from purple to yellow after incubation at 58-60°C, indicating process failure. Absence of color change (medium remains purple) indicates proper sterilization process. To validate the process and ensure sterilization quality, it is recommended to place the biological indicator in a challenge device (DDP) with each load of medical devices or implants.

VALIDITY: 24 months.

DISPOSAL: Dispose of according to the current sanitary regulations in your country.

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