# Durevole – Pino de fibra de vidro

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
# Durevole – Glass fiber posts

## INSTRUCTIONS FOR USE

**DESCRIPTION:** DUREVOLE Glass fiber posts are double intraradicular conical posts. The post has an elasticity module similar to that of dentin, providing a more homogeneous mechanical unit as anchorage for coronary reconstruction, enabling the distribution of stress more evenly along the root. For information on product measurements, see table at the end of the package insert.

**STORAGE:** Store the product in its original, closed packaging, in a clean, dry place with a relative humidity below 95% and at room temperature up to 50°C.

**VALIDITY:** 05 (five) years after the date of manufacture. Single use product. Prohibited reprocessing.

**COMPOSITION:** Glass fibers, epoxy resin and radiopacifying agent.

**INDICATIONS:** Intraradicular posts are indicated to assist in the retention of extensive coronary restoration in endodontically treated teeth.

**CONTRAINDICATIONS:** The post is contraindicated for cases of patients with parafunctional habits, coronary remnants of less than 2 mm, cases of extensive rehabilitation, teeth without adequate endodontic treatment and / or patients with hypersensitivity to the components of the product.

**POST CLEANING:** The cleaning of the post must be carried out with 70% alcohol.

### DIRECTIONS FOR USE:

1. Select the intraradicular post according to the clinical and radiographic characteristics of the intraradicular canal.
2. Perform absolute isolation of the operative field.
3. Using a heated instrument, remove 2/3 of the root canal filling material, preserving 3 to 5 mm of the obturation of the apical third.
4. Prepare the root canal with the use of largo burs.

5. Position the intraradicular post inside the canal and perform a periapical radiography to confirm the adaptation of the post to the root canal.

6. With the use of diamond drill at high speed and under irrigation, perform the cutting of excess coronary post if necessary.

7. Confirm the adjustment of the post inside the intraradicular canal, then clean the post with alcohol and dry.

8. Apply a layer of silane over the entire surface of the post and wait at least 1 (one) minute for complete drying.

9. Etch the root canal for 15 seconds with 37% phosphoric acid (Indusbello) for 15 seconds.

10. Remove the phosphoric acid with plenty of water and dry with an air jet and absorbent paper cones.

11. Apply the adhesive system (according to the manufacturer's instructions) over the entire length of the root canal.

12. Fill the root canal with resin cement and insert the post into the entire prepared canal, keep post it completely covered with cement.


13. Remove excess resinous cement and proceed with photopolymerization, if necessary, according to the cement manufacturer's guidelines.

14. Perform coronary tooth reconstruction.

**WARNINGS:** This product is for dental use only and must be handled in accordance with the instructions for use. Single use product. Prohibited reprocessing. Do not use the product if it is out of date. Observe the storage conditions.

**PRECAUTIONS:** The use of PPE is required to handle the product. The cutting of the post must necessarily be carried out with the use of diamond drills in high rotation and under irrigation. The use of scissors, tweezers and / or pliers is contraindicated for cutting the post. The post must be free of any kind of dirt.

**ADVERSE EFFECTS:** No adverse effects are expected if the product is used in accordance with the instructions.

Durevole Model nº	Colors	A (mm) Height	B (mm) cervical diameter	C (mm) apical diameter
0,3	Green	19	1,25	0,45
0,5	Purple	19	1,40	0,65
1,0	White	19	1,60	0,85
2,0	Yellow	19	1,80	1,05
3,0	Blue	19	2,00	1,25

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