Protetor de Cerdas

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Bristle protector

USE INSTRUCTIONS

INDICATION: For toothbrush bristle protection.

INSTRUCTIONS: Open the bristle protector pressing the "animal's ears" backwards. Place your toothbrush's head at the opening, with bristles tilted out. Then close the Bristle Protector. If you wish, use both silicone suction cups to attach the Protector to the wall.

CLEANING: Use a dry sponge, running water and mild detergent. You can finish cleaning it up by drying the bristles with paper towel or a clean and dry cloth.

STORAGE: Keep the product in closed cabinets protected from dust and insects, and also kept from direct sun and rain.

MATERIAL: ABS.

WARNING: This product is not a toy. Product is intended to be handled by 17+

year-olds and dental professionals. Keep away from children and pets, because it contains small parts that may cause injury.

EXPIRATION DATE: 05 years from manufacturing date, considered a non-violated packaging. Use life of product depends on how it's been handled while in use. Continuous use and sterilization cycles may cause material fatigue and limit its useful life. Verify the condition of this product before each use and purchase another one if you find excessive fatigue and risk of failure.

DISPOSAL OF PRODUCT: Disposal of this product must follow local environmental rules.

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